ENDOVUE™ 24

High-Definition Surgical Visualization





Fully compliant for medical use in surgery, the 24″ EndoVue™ from NDSsi offers a high quality, low-cost solution for surgical imaging environments. Its compact and lightweight design makes it very maneuverable and easy to install on a surgical cart or boom arm. The new 24″ EndoVue™ is an affordable alternative to a full-featured surgical display but can still accommodate high-definition signals from a variety of medical imaging sources, including endoscopes, ultrasound, PACS, and vital signs. It also features a fanless cooling design to minimize the risk of spreading airborne contaminates near the sterile field.

- Purpose-Built for Use in the OR
 - High Quality, Low-cost Solution
 - Lightweight, Highly-Maneuverable
 - Programmable User Preferences
 - Fanless Cooling Design

ENDOVUE™ 24

High-Definition Surgical Visualization



PRODUCT DETAILS

- Large, 24" Diagonal Size
- > HD Imaging Compatible
- High Resolution (1920 x 1200)
- Programmable User and Modality Preferences
- Fanless Cooling, Sealed Front Enclosure
- Multiple Analog & Digital Video Inputs
- Safe Low Voltage 24 VDC Remote Power Supply
- Compact and Lightweight
- Flexible Standard VESA Mounting

SUPPORTED DEVICES/APPLICATIONS

SOURCEAPPLICATIONDVI, HD-SDIHD Endoscopy, PACS, Vital Signs, Room CameraHD-RGBS, HD-YP,P,
RGBS, YP,P,P, SDIHD EndoscopyS-Video, CompositeSD Endoscopy, SD Endoscopy, UltrasoundVGAPACS, Endoscopy, Ultrasound, Angiography

COMPLIANCE & CERTIFICATIONS

UL 60601-1, EN 60601-1, EN 60601-2, MDD 93/42/EEC, CE, CAN/CSA C22.2 NO. 601.1, Class 1 Medical Device, RoHS, FCC Part 15





ABOUT NDS SURGICAL IMAGING

At NDS Surgical Imaging (NDSsi), our vision is built upon the unique ability to supply the medical profession with technologies that enable practitioners to *integrate images* from many sources, *visualize a defect* or condition within the body and to *perform an intervention* that improves the quality of care. Founded in 1996, NDSsi is a privately held company in Silicon Valley with US-based manufacturing, R&D facilities and global sales and service offices. As the company moves into the next decade of innovation, the core focus remains to deliver cutting-edge visualization display products and informatics appliances that enable new clinical imaging solutions and fuel collaborative medicine. All NDSsi products are medical grade devices certified under the most stringent safety and environmental standards.

Learn more about NDS products at: www.ndssi.com



North American Corporate

5750 Hellyer Avenue San Jose, CA 95138 (USA) Tel: 408 776 0085 Toll Free: 866 637 5237 Fax: 408 776 9878 Email: info@ndssi.com

U.S. East Coast Office

Dome Products 400 Fifth Avenue Waltham, MA 02451 United States (USA) Tel: 781 895 1155 Fax: 781 895 1133

European Offices

Nijverheidscentrum 28 2761 JP Zevenhuizen (ZH) The Netherlands Tel: + 31 180 63 43 56 Fax: + 31 180 63 21 91 Fmail: info@ndssi.nl

Asia Pacific Offices

Takanawa Kaneo Bldg., 6F 3-25-22 Takanawa Minato-ku Tokyo, Japan, 108-0074 Tel: + 81 3 5475 1835 Fax: + 81 3 5475 2835 Email: info@ndssi.jp

© 2009 NDS Surgical Imaging. Features and specifications are subject to change without notice. This product is capable of displaying radiology (PACS) images for reference purposes only. Products may not be available in all markets and are subject to the regulatory or medical practices that govern individual markets. Always refer to the package insert, product label and/or user instructions before using any NDSsi product. Contact your NDSsi representative for more information. All trademarks are property of their respective owners.